

## 3-WAY STOPCOCK WITH EXTENSION LINE

**Variants :** 7 cm (Code: 1901) | 10 cm (Code: 1902) | 15 cm (Code: 1903) | 25 cm (Code: 1904) | 50 cm (Code: 1905)  
100 cm (Code: 1905) | 150 cm (Code: 1906)

### Overview

High quality and reliable 3-way stopcock with extension tube for use with any infusion/transfusion devices and blood pressure monitoring. It also provides an alternative channel for introduction of medicine to the patient with ease of use.

### Features

- Individually packed and ready to use
- Disposable, sterile and pyrogenic
- Three arm handle allows for smooth movement to change fluid direction for facilitating smooth and easy operations  
arrow mark indicating direction of flow
- Transparent body - allows easy observation of flow
- Kink-resistant tubing with smooth inner surfaces for laminar, non-turbulent flow.
- Minimum residual fluid volume
- Pressure endurance of > 40 bar

