

## FLEXIBLE CATHETER MOUNT


**Variants** : Without Luer port (Code:4140) | With 4.3 Luer Port (Code:4120) | With 7.6 Luer Port (Code:4130)

### Overview

Catheter mounts are available in corrugated tubing, smooth bore tubing and extendible tubing. A variety of patient elbows are available for gas monitoring, airway suctioning and fibre optic bronchoscopy.

### Features

- Low dead space
- Corrugated tubing and smooth bore tubing
- Available in a wide range of configurations and can be further customised by attaching associated products including filters, face masks etc.



### Specifications

- Connectors : 22 M / 15 F patient connection and 22 F ventilator connection
- Luer port : 4.3 mm / 7.6 mm
- Length : 15 cm