

YANKAUER SUCTION SET

Variants : Standard tip (Code: 1500)

Overview

The suction set can be used for blood and other fluids removal during open surgery while the product is connected to an external source of suction and its suctioned volume counted as blood loss during surgery.

Features

- Individually packed and ready to use
- Disposable, sterile and non-pyrogenic
- Ergonomic transparent handle offers a firm grip
- Handle with 4 side suction ports smooth atraumatic tip to withstand negative pressure and ensure continuous and easy suction of fluid
- Kink resistant transparent tubing makes draining visible
- Connection adaptability to vacuum or suction tubes
- Control vent sealed with removable sleeve
- 2 meter tube

